

Wall-mounted fold-up wheelchair scale

seca
Precision for health

seca XLine Fold

- + Folds up flat to save critical floor space in tight, high-traffic care areas
- + Gas-spring powered deployment eliminates heavy lifting, pinching, or strain – keeping nurses safe from the start
- + Oversized platform with low-incline ramps: 43.3" x 43.3" surface easily accommodates bariatric wheelchairs and transport chairs
- + Built to handle bariatric demands: supports up to 1,100 lbs with rock-solid stability for safe weighing
- + Tare-by-Scanning function instantly subtracts the stored weight of wheelchairs – streamlining workflow and reducing manual steps
- + Recall function, re-weigh function, and multi-side access simplify weighing in fast-paced environments
- + Enhanced ID-Display™ with improved user guidance, smart tare handling and optional EMR functions such as user and patient identification
- + Send weight data directly to the medical record (EMR) – eliminating manual entry and reducing weight-based dosing errors



Image contains barcode scanner from a third-party vendor.

seca XLine Fold

Wall-mounted fold-up wheelchair scale



Unfold Maximum Capability – Without Sacrificing Space

Delivering a full-size, bariatric-ready platform right at the point of care, the seca XLine Fold enables efficient, on-the-spot weighing for patients in wheelchairs, transport chairs, or walkers. Its fold-up design preserves floor space and minimizes the need for patient transfers – saving time, reducing strain, and improving workflow in high-traffic areas.

Model number	620
Capacity	500 kg / 1,100 lbs
Graduation	100 g / 0.2 lb
Dimensions (WxHxD)	1,794 x 1,677 x 1,497 mm / 70.6 x 66 x 58.9 in
Dimensions when folded (WxHxD)	1,794 x 1,677 x 157 mm / 70.6 x 66 x 6.2 in
Platform (WxHxD)	1,100 x 57 x 1,100 mm / 43.3 x 2.2 x 43.3 in
Net weight	115 kg / 331 lbs
Power supply	Power adapter
Functions	Display of user and patient name or ID, Display date of birth patient, TARE, Pre-TARE, HOLD, Auto-HOLD, CLEAR, Auto-CLEAR, RESET, kg/lbs switch-over, kg/lbs switch-over lock, cm/inch switch-over, Adjustable damping, BMI, BSA (Body surface area), Standby mode, Recall
Interfaces	Ethernet, Wi-Fi, RS-232 (optional)
Medical Device Regulation (MDR) certified	



Effortless Operation with Built-In Efficiency

Secure wall-mounting and gas-spring powered deployment eliminates heavy lifting, pinching, and strain – keeping nurses safe from the start. The ultra-low profile and third ramp enable smooth, multi-directional access for all patient types. Smart features include a Recall function to retrieve the last recorded weight and Pre-TARE memory to store up to three preset deductions – ideal for commonly used mobility aids.



Seamless Integration for Ultimate Safety

Validated for integration with leading EMR systems like Cerner, Epic, MEDITECH, and HL7, the seca XLine Fold eliminates manual weight entry – reducing documentation errors and safeguarding weight-based medication dosing. The intelligent ID-Display™ supports nurse-to-patient identification by showing key information such as the patient's name, user ID, and date of birth directly on the scale.

seca headquarters
(responsible for CIS & Africa)
+49 40 20 00 00 0
info@seca.com
seca gmbh & co. kg
Hamburg, Germany

seca middle east
+971 50 410 2319
info.ae@seca.com
Dubai
United Arab Emirates

seca asia pacific
+60 12 2170554
info.aspac@seca.com
Kuala Lumpur
Malaysia

seca.com/XLine